



Flamingo Instructions for Use

Table of Contents

Introduction	3
Product overview	4
Safety precautions & troubleshooting	5
Assembling MCU and Camera module	6
Flamingo MCU indicator	7
Buttery charging	8
Download Flamingo software	9
Technical specifications	10
Warranty	11

Introduction

Thank you for choosing **adme**tec

Intended Use:

Flamingo is an integrated HD camera and light system designed for use in dental and medical practices. It attaches to magnification loupes to assist practitioners in documenting and illuminating procedures. This innovative device enables real-time display and recording of detailed images, enhancing both patient care and educational opportunities for medical professionals.

About Flamingo

• Lightweight (just 19 grams!)

Designed to be a lightweight ergonomic extension of your loupes.

• Integrated light and camera

Proper lighting is critical both for focus and quality image capture.

• Live broadcast via Wi-Fi to the computer

Real time display of what you are doing, while you do it (delay of 0.05 seconds).

• Full HD video

Ultra-smooth, high-quality 1080p | 30fps, vivid, "living color" video recording.

• Unrestricted movement

With Admetec's Main Control Unit (MCU) on your belt or in your pocket. For unparalleled convenience and user-friendliness.

• Admetec proprietary software

Enables video editing including taking quality screen captures for real-time collaboration and documentation.

Product overview



Flamingo Insert Set Please note that the image is for illustration.

Your Flamingo comes with:

- 1. WiFi USB dongle (for Windows)
- 2. Camera & light module
- 3. Camera module cable to MCU
- 4. Charging stand
- 5. MCU with 2 batteries
- 6. AC/DC adaptor
- 7. Wireless Pedal
- 8. Power cable to the AC/DC adaptor
- 9. Pedal Type-C charging cable
- 10.Screwdriver
- 11. Frame-temple adaptor for cable
- 12. Clothing cable-clip
- 13. Wi-Fi USB dongle (for Mac iOS)

*The WiFi dongle will be supplied according to user's order (Windows or iOS)

Safety Precautions

Do not use any third-party AC/DC adaptor with Flamingo battery charger stand, even if they look identical. Third-party chargers may not be compatible with the battery charger and could damage it.

Prevent the MCU from overheating. Avoid direct sunlight or heated surfaces.

Before using the device, verify that it is assembled and connected as described on page 5.

When used with non-Admetec loupes, use the appropriate adaptor to ensure maximum focus and stability when using Flamingo camera.

Do not use in proximity to explosive or highly flammable materials such as: alcohol, methanol, acetone, etc.

Please note that any serious incident should be reported directly to the manufacturer with relevant documentation describing the event.

Troubleshooting

Issue	Troubleshooting
The light does not work steadily	Check cable. Contact your local
	dealer for instructions
The video does not record or show video	Contact local dealer or
image	flamingo@admetec.com
Low brightness	Use the brightness adjustment
	buttons [Page 6]
Battery overheats regularly	Contact local dealer
Buttons are stuck	Contact local dealer
Flamingo application issues email:	flamingo@admetec.com

Assembling MCU and Camera Module

To insert the battery into the MCU socket.

- 1. Remove protective foil from the battery contacts. [Fig. 01]
- 2. Insert the battery into the MCU socket. [Fig. 02]



Fig. 01 – Flamingo battery



Connect the MCU and camera module with the provided cable.

- 1. Plug the USB Type C side into the MCU port [Fig. 03]
- 2. Plug the Micro USB cable type into the camera module port. [Fig. 04]



Fig. 03 – Type C cable insert into MCU



Fig. 04 – Camera module with Micro USB

Caution



Do not use any third-party AC/DC adaptor with Flamingo battery charger stand. Third-party chargers damage the battery charger stand.

Caution



Prevent MCU from overheating.

To check battery status, press the "Check" button (see Page 06).

Flamingo MCU Indicators

Each Flamingo kit includes one MCU and two batteries which serve as a power source and control panel for the camera and light module.



Item/Button	Description	Action
CAMERA	Camera Power	To turn ON, press the 'ON/OFF' button.
ON/OFF PC REC	ON/OFF button	To turn OFF, press & hold for 3 sec' the
		'ON/OFF' button.
CAMERA	REC Button	Press the 'REC' button to enter recording
ONIOFF PC REC		mode and recording.
CAMERA	LED Indicators	Blinking blue = MCU is loading
ON/OFF PC REC		Steady blue = Device is ready
60		Green led = Linked to PC/Laptop
		Blinking red = Recording mode ready
		Constant red = Recording live mode
LIGHT	Light Power	To turn ON or OFF, press the 'ON/OFF'
+ ON/OFF	ON/OFF button	button.
LIGHT	Light Intensity	Press the '+' or '-' to increase or
+ ON/OFF	controls	decrease the light intensity.
	Battery status	To check battery level press & hold the
LOW FULL		'Check' button.

▲	Light Power	To turn ON or OFF, press the 'ON/OFF'
*	ON/OFF cable	button.
	switch	

Battery charging

 Connect the AC/DC adaptor output to the Flamingo charging stand. Make sure the power cable is connected to the AC/DC adaptor and plugged to the wall power outlet.



2. Insert the battery to the charging stand socket as shown below.



Average charging time for an empty battery is 3.0 hours.

Charging stand - Light Indicators

Blue color indicator	Power connected
Red color indicator	Battery is charging
Green color indicator	Battery is fully charged

Downloading the software

To set up your Flamingo you are required to download Admetec's proprietary software on your device/PC and follow the instructions in the user manual on the website.

Download the Flamingo software and the User manual from the following link: <u>https://www.admetec.com/flamingo</u>

Minimum PC Requirements:

- OS: Windows 10 or above
- CPU: Intel Core i5 or better or AMD Ryzen 5 or better
- GPU: An integrated graphics card
- RAM: Minimum of 8GB
- Storage: Minimum of 250GB SSD
- USB: one free USB 2.0/3.0 port

Technical Specifications

Battery chargerInput voltage100 - 240 V ACOutput voltage5 V DC

Main Control Unit parameters

Internal battery voltage:	3,7 V
Battery capacity:	6.5 Ah
Full battery charge time	2.5 hrs
Average working time (with camera & light together)	> 4 hrs
Average working time (with camera or light only)	> 8 hrs
Light digital dimmer	32-step
Height:	101 mm
Width:	80 mm
Depth: (without belt clip)	28 mm
Weight: (with battery)	229 gr

Camera & light source

Camera lens	25mm IR 1/2" 5MP
Illuminance	70,000 lx
Color temperature	5750 °K
Light diameter	17 mm
Weight (without cable)	19 gr
Spotlight size diameter at 350mm	90mm
Resolution	1080p FHD
Sensor size	1/2.7 inch

Flamingo pedal

Voltage	3,7 V
Size	W: 52 mm / D: 98 mm / H: 25 mm
Weight	63 gr
вт	Bluetooth Model 4.2/5
Charger Port	USB Type C

Warranty

Admetec's warranty obligations for the Flamingo are limited to the terms set forth below:

Admetec warrants the Flamingo product against defects in materials and workmanship for a period of two (2) years for the Flamingo set and one (1) year for all other supplied parts from the date of original purchase ("Warranty Period").

If a defect arises and a valid claim is received by Admetec within the Warranty Period, Admetec will (1) repair the product at no charge, using new or refurbished replacement parts, (2) exchange the product with a product that is new, or which has been manufactured from new or serviceable used parts and is at least functionally equivalent to the original product.

If a defect arises and a valid claim is received by Admetec after the first ninety (90) days of the Warranty Period, a shipping and handling charge will apply to any repair or exchange of the product undertaken by Admetec.

When a product or part thereof is exchanged, any replacement item becomes your property, and the replaced item becomes Admetec's property.

Exclusions and Limitations

This warranty does not apply to: (a) to damage caused by accident, abuse, misuse, misapplication, or non-Admetec products; (b) to damage caused by service performed by anyone other than Admetec; (c) to a product or a part that has been modified without the written permission of Admetec; or (d) if any Admetec warranty label has been removed or defaced.

CE	CE Mark indicating that the device complies with Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices
	Admetec Solutions Ltd. Yokhanan HaSandlar St. 18 Haifa, 3296218, Israel
EC REP	Obelis Group Bd Général Wahis, 53, Brussels, B-1030 Belgium
MD	Medical Device
$\underline{\land}$	Caution! Consult accompanying documents
LOT	Batch number
REF	Catalog number



If you have purchased Admetec-branded electrical or electronic products in the EU and are intending to discard these products at the end of their useful life, please do not dispose of them with your other household or municipal waste. Admetec has labeled its branded electronic products with the WEEE Symbol to alert our customers that products bearing this label should not be disposed

of in a landfill or with municipal or household waste in the EU. Instead, contact the local e-waste recycling center.

IM4000005 Rev.A

Do you need customer service?

Please contact your local distributor.



www.admetec.com