

The art of seeing
what others can't.



Handcrafted Loupes, Lights & Camera

admetec
We've got your back





admetec ergo

Your body wasn't designed
to bend over patients
for hours on end.

Ergonomic loupes

Admetec Ergo comes in 6 magnifications.

Designed specifically for your working distance, to the centimeter.



3.0x

Working distance range
35-70 cm

Weight
51gr



4.0x

Working distance range
35-70 cm

Weight
52gr



5.0x

Working distance range
35-70 cm

Weight
53gr



6.0x

Working distance range
35-70 cm

Weight
53gr



7.5x

Working distance range
40-70 cm

Weight
54gr

The highest ergonomic magnification power available
for the pinpoint accuracy you need.

10x magnification power

Working distance range
40-70 cm

Weight
58gr

Ergo prism technology
developed by Admetec

Lightweight, Polycarbonate
lenses for durability
and a premium image







Ocular lenses from the
Japanese world leader in glass
lenses manufacturing, FUJITOK®

What style makes your heart sing?

You're still YOU while you are at work. That's why we made sure to provide a wide variety of frames that allow you to express your unique personality and work in comfort.



Blues

-  Navy Blue
-  Bronze & Light Blue
-  Black
-  Pink







Indie

-  Champagne
-  Orange Lava
-  Midnight Blue
-  Black






Soul

-  Pink
-  Brown & Light Blue
-  Burgundy & Gold
-  Black



Jazz

-  Red & Black
-  Orange & Black
-  Blue & Black



Morriz of Sweden®

-  Red
-  Green
-  Blue
-  Black



Progear®

-  Red & Metallic Brown
-  Deep Green & White
-  White & Blue
-  Black



Bollé®

-  Gun Brown





We've got **your back.**

Admetec designs, develops, manufactures, and distributes handcrafted, custom-fitted magnifying loupes and LED headlights as well as an integrated camera and headlight. We specialize in creating products that are exceptionally lightweight, convenient to use, and ergonomic – enabling you to focus on being effective while we focus on your comfort.

MM4000001



+972-4-989-0747 info@admetec.com www.admetec.com